

Respiratory Therapy

Disposable Resuscitator

Description



Adult

SW71100B Disposable Resuscitator Adult

Capacity of resuscitator	1600ml \pm 10%
Material of resuscitator	PVC
Pop-off valve	60cm H ₂ O
Disposable reservoir bag	1600ml
Disposable mask	#5
Oxygen tube	2m
Air Way	80mm-100mm



Child

SW71200A Disposable Resuscitator Child

Capacity of resuscitator	500ml \pm 10%
Material of resuscitator	PVC
Pop-off valve	40cm H ₂ O
Disposable reservoir bag	1600ml
Disposable mask	#3
Oxygen tube	2m
Air Way	60mm-80mm



Infant

SW71300A Disposable Resuscitator Infant

Capacity of resuscitator	280ml \pm 10%
Material of resuscitator	PVC
Pop-off valve	40cm H ₂ O
Disposable reservoir bag	1000ml
Disposable mask	#1
Oxygen tube	2m
Air Way	40mm-60mm

Respiratory Therapy

Reusable Resuscitator



Adult

SW72100A Reusable Resuscitator Adult

Capacity of resuscitator	1600ml ± 10%
Material of resuscitator	PC, Silicon
Pop-off valve	60cm H ₂ O
Disposable reservoir bag	2500ml
Reusable mask	#3
Oxygen tube	2m
Air Way	80mm-100mm



Child

SW72200A Reusable Resuscitator Child

Capacity of resuscitator	500ml ± 10%
Material of resuscitator	PC, Silicon
Pop-off valve	40cm H ₂ O
Disposable reservoir bag	2500ml
Reusable mask	#2
Oxygen tube	2m
Air Way	60mm-80mm



Infant

SW72300A Reusable Resuscitator Infant

Capacity of resuscitator	280ml ± 10%
Material of resuscitator	PC, Silicon
Pop-off valve	40cm H ₂ O
Disposable reservoir bag	1000ml
Reusable mask	#0
Oxygen tube	2m
Air Way	40mm-60mm