













Electrosurgical Pencil

Reusable Pencil

Ø2.38mm

Description	Package	Ref. No.
Reusable hand control electrosurgical pencil, button switch, 70mm non-stick blade, 3m silicon cable, 3-pin	 =80	SW11200-ASBCS
		
Reusable hand control electrosurgical pencil, button switch, 70mm non-stick blade, 5m silicon cable, 3-pin	 =80	SW11200-ASBCS-BE501
		
Reusable hand control electrosurgical pencil, button switch, 70mm non-stick blade, 3m silicon cable, 4mm one pin	 =80	SW11200-DSBCS
		
Reusable hand control electrosurgical pencil, button switch, 70mm non-stick blade, 3m silicon cable, 3-pin, *Fixed cap	 =80	SW11200-T03
		

Button Switch

70mm non-stick blade

Ref. no.	Cable	Connection	Package	Carton	Accessories
SW11200-ASBCS	3m	3-pin	Single	80pcs	
SW11200-ASBCS-BE501	5m	3-pin	Single	80pcs	
SW11200-DSBCS	3m	4mm erbe	Single	80pcs	
SW11200-T03	3m	3-pin	Single	80pcs	*Fixed cap

70mm stainless blade

Ref. no.	Cable	Connection	Package	Carton	Accessories
SW11200-ASBCS-B	5m	3-pin	Single	80pcs	

Electrosurgical Pencil

Reusable Pencil

Ø2.38mm

Description

Package

Ref. No.

Reusable foot control electrosurgical pencil, 70mm non-stick blade, 3m silicon cable, 8mm adapter



=80

SW11100

Type: SW11100-CSFAS



Foot Control

70mm non-stick blade

Ref. no.	Cable	Connection	Package	Carton	Accessories
SW11100	3m	4mm pin	Single	60pc	4/8mm adapter
SW11100-CSFAS	3m	8mm adapter	Single	60pc	