


Smoke Evacuation Pencil

ShinEvac[®] AIO (All-In-One)

Description	Package	Model No.
Disposable hand control smoke evacuation pencil, telescopic tube, with 360° patented swivel design & double package with EO sterilization. Accessories included a clean-edge holster and connectors.	 =20	SW12200-APBCSB36
		SW12200-APBCSB56

360° patented swivel design.



510K NO. : K213317

Accessories

Clean-Edge Holster



1.8**1.2**5.3"
(45*30*135mm)

Connectors



0.4"-0.4"
(10mm-10mm)
&
0.4"-0.9"
(10mm- 22mm)

Ordering Information

Model no.	Cable Length	Insulated Electrode	Smoke Suction Tube
SW12200-APBCSB36	9.8' (3M)	Active 0.6"/2.6"(15/70mm) Non-Stick blade	Extendable 0.98"-5.9" (25mm-150mm)
SW12200-APBCSB56	16.4' (5M)		




升望科技股份有限公司 Shining World Health Care Co., Ltd.


Tel: +886-2-22900966 Fax: +886-2-22903966 www.shinmed.com export@shinmed.com

Smoke Evacuation Pencil

ShinEvac® Classic

Description	Package	Model No.
Disposable hand control smoke evacuation pencil, 2 stage adjustable extension tube, with 360°patented swivel design & double package with EO sterilization. Accessories included a clean-edge holster and connectors.	 =20	SW12200-APBCSE36
		SW12200-APBCSE36G
		SW12200-APBCSE56
		SW12200-APBCSE56G

360° patented swivel design.



510K NO. : K213317



台灣精品
TAIWAN EXCELLENCE
2019

Accessories

Clean-Edge Holster



1.8**1.2**5.3"
(45*30*135mm)

Connectors



0.4"-0.4"
(10mm-10mm)
&
0.4"-0.9"
(10mm- 22mm)

Ordering Information

Model no.	Cable Length	Insulated Electrode	Smoke Suction Tube
SW12200-APBCSE36	9.8'	Active 0.6"/2.6" (15/65mm) Non-Stick blade	1.2" (30mm)
SW12200-APBCSE36G (Non-PVC)	(3M)		
SW12200-APBCSE56	16.4'		
SW12200-APBCSE56G (Non-PVC)	(5M)		