


Smoke Evacuation Pencil

ShinEvac® AIO (All-In-One)

Description	Package	Ref. No.
Disposable hand control smoke evacuation pencil, telescopic tube, active 15mm/70mm non-stick blade electrode, 3M & 5M cable, 10-10mm & 10-22mm connector, with 360° patented swivel design.	 =20	SW12200-APBCSD01

360° patented swivel design.



衛署醫器製字第 003801 號

CE 2460

Ordering Information

Ref. no.	Cable Length	Electrode	Smoke Suction Tube
SW12200-APBCSB02	3M	Active 15/70mm Non-Stick blade	Extendable (25mm-150mm)
SW12200-APBCSD01		Active 15/70mm stainless blade	
SW12200-APBCSA01	5M	Active 15/70mm Non-Stick blade	
SW12200-APBCS		Active 15/70mm stainless blade	

Accessories

Electrodes



Supports
70-180mm
electrodes

Connectors





10mm-10mm

10mm-22mm

Smoke Evacuation Pencil

ShinEvac® Classic

Description	Package	Ref. No.
Disposable hand control smoke evacuation pencil, 2 stage adjustable extension tube, active 15mm/70mm non-stick blade electrode, 3M cable, 10-10mm & 10-22mm connector, with 360° patented swivel design.	 =20	SW12200-APBCSE04



360° patented swivel design.

衛署醫器製字第 003801 號


CE 2460



Ordering Information

Ref no.	Cable Length	Electrode	Smoke Suction Tube
SW12200-APBCSE04	3M	Active 15/65mm Non-Stick blade	30mm
SW12200-APBCSE03		Active 15/65mm stainless blade	
SW12200-APBCSE91 (Non-PVC)		Active 15/65mm Non-Stick blade	
SW12200-APBCSE08	5M	Active 15/65mm Non-Stick blade	80mm
SW12200-APBCSE07		Active 15/65mm stainless blade	
SW131201117S-L		Active 15/120mm Non-Stick blade	


Electrodes



Supports
50mm-70mm &
120mm
electrodes

10mm-10mm

Accessories Connectors



10mm-22mm

X-Ray Detectable Extension Tube

80mm 

30mm 