



















Electrosurgical Accessories

Bipolar Forceps Cable

<p>Male Plug (Generator)</p> 		
	<p>Female Plug (Instrument)</p> <p>ø 4mm x 2 prong pin</p>	<p>28.8mm</p> <p>Valleylab/Bowa/Aesculap</p> <p>Lamidea/Olympus/Soering</p>
 <p>European type</p>  <p>3.8mm</p> <p>6.7mm</p> <p>European Type 3.8mm Flat Plug Length: 3M</p>	<p>SW4643002A(Grey)</p>	<p>SW4643001E(Grey)</p>
 <p>European type</p>  <p>4.2mm</p> <p>6.7mm</p> <p>European Type 4.2mm Flat Plug Length: 3M</p>	<p>SW4643004A(Blue)</p> <p>SW4643004B(Blue)</p>	<p>SW4643004D(Blue)</p> <p>SW4643004G(Blue)</p> <p>SW4653001A(Blue)</p> <p>SW4653001B(Blue)</p>
 <p>US type</p>  <p>5.50mm</p> <p>2</p> <p>US Type ø 2.0mm Angled Plug Length: 3M</p>	<p>SW4613005A(Blue)</p> <p>SW4613005E(Blue)</p> <p>SW4613003E(Grey)</p>	<p>SW4613005B(Blue)</p> <p>SW4613005D(Blue)</p> <p>SW4613003F(Grey)</p>
  <p>US Type ø 2.0mm Booted Plug Length: 3M</p>	<p>SW4613003I(Grey)</p>	<p>SW4613003H(Grey)</p>
  <p>Aesculap Type Booted Plug Length: 3M</p>	<p>SW4643007B(Blue)</p> <p>SW4643007C(Grey)</p>	

Bipolar Forceps Cable

<p>Male Plug (Generator)</p>	
	<p>ø 4mm x 2 prong pin</p>
<p>Female Plug (Instrument)</p>  <p>US Type ø 2mm Non Booted plug Length: 3M</p>	<p>SW4613003I-01(Grey)</p>
  <p>Aesculap Type Non Booted Plug Length: 3M</p>	<p>SW4643007D(Grey)</p>
 <p>2.0mm x 2 Flat Plug</p>	<p>SW4613004B(Grey)</p>

Description

Item Number	Feature	Color	Type	Description		Length	Material	Packaging	Sterile
				Male Plug	Female Plug				
				(Generator side)	(Instrument side)				
SW4613005A-7	Disposable	Blue	US	4mm x 2 prong pin	2.0mm flat plug	3 M	MPVC	Single packaging, 100 pcs / ctn	Y
SW4613005B	Disposable	Blue	US	28.8mm 2 prong pin	2.0mm flat plug	3 M	MPVC	Single packaging, 100 pcs / ctn	Y
SW4613005D	Disposable	Blue	US	28.8mm 2 prong pin	2.0mm flat plug	3 M	MPVC	Bulk, 250 pcs / ctn	N
SW4613005E	Disposable	Blue	US	4mm x 2 prong pin	2.0mm flat plug	3 M	MPVC	Bulk, 250 pcs / ctn	N
SW4643004A	Disposable	Blue	European	4mm x 2 prong pin	4.2mm flat plug	3 M	MPVC	Single packaging, 100 pcs / ctn	Y
SW4643004B	Disposable	Blue	European	4mm x 2 prong pin	4.2mm flat plug	3 M	MPVC	Bulk, 250 pcs / ctn	N
SW4643004D	Disposable	Blue	European	28.8mm 2 prong pin	4.2mm flat plug	3 M	MPVC	Bulk, 250 pcs / ctn	N
SW4643004G	Disposable	Blue	European	28.8mm 2 prong pin	4.2mm flat plug	3 M	MPVC	Single packaging, 100 pcs / ctn	Y
SW4653001A	Disposable	Blue	European	2.0mm MARTIN plug	4.2mm flat plug	3 M	MPVC	Single packaging, 100 pcs / ctn	Y
SW4653001B	Disposable	Blue	European	2.0mm MARTIN plug	4.2mm flat plug	3 M	MPVC	Bulk, 250 pcs / ctn	N
SW4613003E	Reusable	Grey	US	4mm x 2 prong pin	2.0mm flat plug	3 M	TPR plug + Silicon cord	Single packaging, 80 pcs / ctn	Y
SW4613003F	Reusable	Grey	US	28.8mm 2 prong pin	2.0mm flat plug	3 M	TPR plug + Silicon cord	Single packaging, 80 pcs / ctn	N
SW4613003H	Reusable	Grey	US	28.8mm 2 prong pin	2.0mm booted plug	3 M	TPR plug + Silicon cord	Single packaging, 80 pcs / ctn	N
SW4613003I	Reusable	Grey	US	4mm x 2 prong pin	2.0mm booted plug	3 M	TPR plug + Silicon cord	Single packaging, 80 pcs / ctn	N
SW4613003I-01	Reusable	Grey	US	4mm x 2 prong pin	2.0 mm non booted plug	3M	TPR plug + Silicon cord	Single packaging, 80 pcs / ctn	N
SW4613004A	Reusable	Blue		4mm x 2 prong pin	2.4mm x 2 flat plug	3 M	TPR plug + Silicon cord	Single packaging, 80 pcs / ctn	N
SW4613004B	Reusable	Grey		4mm x 2 prong pin	2.0mm x 2 flat plug	3 M	TPR plug + Silicon cord	Single packaging, 80 pcs / ctn	N
SW4643001E	Reusable	Grey	European	28.8mm 2 prong pin	3.8mm flat plug	3 M	TPR plug + Silicon cord	Single packaging, 80 pcs / ctn	N
SW4643002A	Reusable	Grey	European	4mm x 2 prong pin	3.8mm flat plug	3 M	TPR plug + Silicon cord	Single packaging, 80 pcs / ctn	N
SW4643007A	Reusable	Grey	Aesculap	28.8mm 2 prong pin	Flat plug	3 M	TPR plug + Silicon cord	Single packaging, 80 pcs / ctn	N
SW4643007B	Reusable	Blue	Aesculap	4mm x 2 prong pin	Flat plug	3 M	TPR plug + Silicon cord	Single packaging, 80 pcs / ctn	N

SW4643007C	Reusable	Grey	Aesculap	4mm x 2 prong pin	Flat booted plug	3 M	TPR plug + Silicon cord	Single packaging, 80 pcs / ctn	N
SW4643007D	Reusable	Grey	Aesculap	4mm x 2 prong pin	Flat non booted plug	3 M	TPR plug + Silicon cord	Single packaging, 80 pcs / ctn	N

