











Electrosurgical Accessories

Bipolar Forceps Cable

Male Plug (Generator)				
Female Plug (Instrument)	ø 4mm x 2 prong pin	28.8mm Valleylab/Bowa/Aesculap Lamideay/Olympus/Soering	24mm Alcon Type	2.0mm Martin/Aesculap Wolf
 <p>European type 3.8mm Flat Plug Length: 3M</p>	SW4643002A(Grey)	SW4643001E(Grey)		
 <p>European type 4.2mm Flat Plug Length: 3M</p>	SW4643004A(Blue) SW4643004B(Blue)	SW4643004D(Blue) SW4643004G(Blue)		SW4653001A(Blue) SW4653001B(Blue)
 <p>US type US Type ø 1.8mm Angled Plug Length: 3M</p>	SW4613003A(Blue)	SW4613003D(Blue)		
 <p>US type US Type ø 2.0mm Angled Plug Length: 3M</p>	SW4613005A(Blue) SW4613005E(Blue) SW4613003E(Grey)	SW4613005B(Blue) SW4613005D(Blue) SW4613003F(Grey)	SW4613005C(Blue)	
 <p>US Type ø 2.0mm Protect Plug Length: 3M</p>	SW4613003I(Grey)	SW4613003H(Grey)		
 <p>Aesculap Type Protect Plug Length: 3M</p>	SW4643007B(Blue) SW4643007C(Grey)	SW4643007A(Grey)		

Description

Item Number	Feature	Color	Type	Description		Length	Material	Packaging	Sterile
				Male Plug	Female Plug				
				(Generator side)	(Instrument side)				
SW4613005A-7	Disposable	Blue	US	4mm x 2 prong pin	2.0mm flat plug	3 M	MPVC	Single packaging, 100 pcs / ctn	Y
SW4613005B	Disposable	Blue	US	28.8mm 2 prong pin	2.0mm flat plug	3 M	MPVC	Single packaging, 100 pcs / ctn	Y
SW4613005C	Disposable	Blue	Alcon	24mm 2 prong pin (Alcon)	2.0mm flat plug	3 M	MPVC	Single packaging, 100 pcs / ctn	Y
SW4613005D	Disposable	Blue	US	28.8mm 2 prong pin	2.0mm flat plug	3 M	MPVC	Bulk, 250 pcs / ctn	N
SW4613005E	Disposable	Blue	US	4mm x 2 prong pin	2.0mm flat plug	3 M	MPVC	Bulk, 250 pcs / ctn	N
SW4643004A	Disposable	Blue	European	4mm x 2 prong pin	4.2mm flat plug	3 M	MPVC	Single packaging, 100 pcs / ctn	Y
SW4643004B	Disposable	Blue	European	4mm x 2 prong pin	4.2mm flat plug	3 M	MPVC	Bulk, 250 pcs / ctn	N
SW4643004D	Disposable	Blue	European	28.8mm 2 prong pin	4.2mm flat plug	3 M	MPVC	Bulk, 250 pcs / ctn	N
SW4643004G	Disposable	Blue	European	28.8mm 2 prong pin	4.2mm flat plug	3 M	MPVC	Single packaging, 100 pcs / ctn	Y
SW4653001A	Disposable	Blue	European	2.0mm MARTIN plug	4.2mm flat plug	3 M	MPVC	Single packaging, 100 pcs / ctn	Y
SW4653001B	Disposable	Blue	European	2.0mm MARTIN plug	4.2mm flat plug	3 M	MPVC	Bulk, 250 pcs / ctn	N
SW4613003D	Reusable	Blue	US	28.8mm 2 prong pin	1.8mm angled plug	3 M	TPR plug + Silicon cord	Single packaging, 80 pcs / ctn	N
SW4613003E	Reusable	Grey	US	4mm x 2 prong pin	2.0mm flat plug	3 M	TPR plug + Silicon cord	Single packaging, 80 pcs / ctn	Y
SW4613003F	Reusable	Grey	US	28.8mm 2 prong pin	2.0mm flat plug	3 M	TPR plug + Silicon cord	Single packaging, 80 pcs / ctn	N
SW4613003H	Reusable	Grey	US	28.8mm 2 prong pin	2.0mm protect plug	3 M	TPR plug + Silicon cord	Single packaging, 80 pcs / ctn	N
SW4613003I	Reusable	Grey	US	4mm x 2 prong pin	2.0mm protect plug	3 M	TPR plug + Silicon cord	Single packaging, 80 pcs / ctn	N
SW4613004A	Reusable	Blue		4mm x 2 prong pin	2.4mm x 2 flat plug	3 M	TPR plug + Silicon cord	Single packaging, 80 pcs / ctn	N
SW4613004B	Reusable	Grey		4mm x 2 prong pin	2.0mm x 2 flat plug	3 M	TPR plug + Silicon cord	Single packaging, 80 pcs / ctn	N
SW4643001E	Reusable	Grey	European	28.8mm 2 prong pin	3.8mm flat plug	3 M	TPR plug + Silicon cord	Single packaging, 80 pcs / ctn	N
SW4643002A	Reusable	Grey	European	4mm x 2 prong pin	3.8mm flat plug	3 M	TPR plug + Silicon cord	Single packaging, 80 pcs / ctn	N
SW4643007A	Reusable	Grey	Aesculap	28.8mm 2 prong pin	Flat plug	3 M	TPR plug + Silicon cord	Single packaging, 80 pcs / ctn	N

SW4643007B	Reusable	Blue	Aesculap	4mm x 2 prong pin	Flat plug	3 M	TPR plug + Silicon cord	Single packaging, 80 pcs / ctn	N
SW4643007C	Reusable	Grey	Aesculap	4mm x 2 prong pin	Flat plug	3 M	TPR plug + Silicon cord	Single packaging, 80 pcs / ctn	N

