

# Ultrasonic Nebulizer

**SW966**

**Home Care**

Description



SW96610	115V / 60Hz
SW96620	230V / 50 HZ
Particle Size	0.5 – 6.0 $\mu$ M
Oscillation Frequency	1.63 MHz
Medication Cup	30 c.c.
Water Tank	180 c.c.
Nebulizing Rate	$\approx$ 1c.c./min
Volume Switcher	Yes
Power consumption	25 Watts
Operating Temperature	10°C to 40°C ( 50°F to 104°F )
Storage Temperature	5°C to 80°C ( 41°F to 176°F )
Unit Weight	1.25 kg
Size	18cm(L)*10.5cm(W)*21cm(H) ( 7.08’’*4.13’’*8.26’’ )
Accessories	Medication Cups * 5 Mouthpiece * 1 Flexible Connection Hose* 2 Adapter * 2 Child Mask * 1 Adult Mask * 1