

Ultrasonic Nebulizer

SW918

Hospital

Description



SW91810	115V / 60Hz
SW91820	230V / 50 HZ
Particle Size	1.0 – 5.0 μ M
Oscillation Frequency	1.7 MHz
Medication Cup	50 c.c.
Water Tank	250 c.c.
Nebulizing Rate	\approx 4 c.c./min
Volume Switcher	Yes
Power consumption	60 Watts
Unit Weight	2.1 kg
Size	30cm(L)*18cm(W)*27cm(H) (11.8''*7.1''*10.6'')
Accessories	Medication Cups * 5 Mouthpiece * 2 75cm Hose* 2 Connector * 4 Child Mask * 1 Adult Mask * 1